

Padding and Lining Luxofoam

GENERAL PRODUCT INFORMATION

Luxofoam is a self-adhesive soft PVC foam with closed cells. It is used as a softening material on orthoses, external immobilization devices and rehabilitation aids.

PRODUCT RANGE

Luxofoam is available in sheets of 50 cm (width) x 97 cm (length) x 3 mm (thickness).

Article Number	Description 1 BC
35110	Luxofoam - Non Perforated
35112	Luxofoam - Perforated

PRECAUTIONS BEFORE USE

- Cut the required pattern before removing the paper from the back.
- Stick Luxofoam on the orthosis either before activating the orthotic fabrication material or after making the orthosis. In the latter case, the extra space which is taken by Luxofoam should be taken into account.
- Luxofoam can still be easily removed immediately after application. After a while, this will be difficult, unless the orthosis is reheated.

MAINTENANCE

- The material must be cleaned daily. Use lukewarm water and disinfecting soap or pre-moistened isopropanol wipes. Rinse well and dry thoroughly.
- Sterilization of the product in an autoclave is impossible. Sterilization by means of gas treatment is possible. Disinfection is possible with alcohol, quaternary ammonium or a solution of commercial disinfection soaps (HAC®, Sterilium®, etc.).
- **Never use solvents.**

WASTE MANAGEMENT

- After use, the product can be disposed of with normal household waste without harming the environment.

ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets, and regulatory information, please visit our website www.orfit.com.

Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries.
This product is a registered trademark from ORFIT INDUSTRIES BV.



ORFIT INDUSTRIES BV

Vosveld 9A | B-2110 Wijnegem | Belgium

T (+32) (0)3 326 20 26

welcome@orfit.com

www.orfit.com



Ref. Nr.	31300
Version	12
Last updated	1/12/2024
Revision date	1/12/2026