

Orfit® Classic

Gauntlet thumb post orthosis

GENERAL PRODUCT INFORMATION

This precut is made out of a low-temperature thermoplastic material and can be molded directly on the patient.

PRODUCT RANGE

This precut is available in different sizes.

Article Number	Color	Thickness (mm)	Size	Perforation Type
35810	Beige	1,6	S	Micro Perforated
35811	Beige	1,6	M	Micro Perforated
35812	Beige	1,6	L	Micro Perforated

INDICATIONS


- Acute MCP arthristis
- Basal joint osteo arthritis
- Rheumatoid arthritis
- Thumb sprain
- Median nerve palsy

CONTRAINDICATIONS

This pattern is NOT SUITABLE for the following applications:

- Fractures of the scaphoide
- Acute De Quervain tendinitis

SIZING



SIZING


Measure the total hand width at the MCP joints of the 2nd to 5th finger.

60-70 mm (2 ¹/₅" - 2 ³/₄") = **SMALL**
 80-90 mm (3 ¹/₈" - 3 ¹/₂") = **MEDIUM**
 100-110 mm (4" - 4 ³/₁₀") = **LARGE**

GENERAL INFORMATION

To take precautionary measures, learn about activation methods , maintenance and waste management read the general Instructions For Use of Orfit Classic.

APPLICATION METHOD



PRECAUTIONS BEFORE USE

ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets, and regulatory information, please visit our website www.orfit.com.

Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries.
This product is a registered trademark from ORFIT INDUSTRIES BV.



ORFIT INDUSTRIES BV
Vosveld 9A | B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

www.orfit.com



Ref. Nr.	35810-35811-35812 ORFIT CLASSIC
Version	27
Last updated	1/12/2024
Revision date	1/12/2026