

## Orfit® Classic

### Rheumatoid arthritis resting orthosis

#### GENERAL PRODUCT INFORMATION

This precut is made out of a low-temperature thermoplastic material and can be molded directly on the patient.

#### PRODUCT RANGE

This precut is available in different sizes.

Article Number	Color	Thickness (mm)	Size	Perforation Type
35840	Beige	3,2	S	Mini Perforated
35841	Beige	3,2	M	Mini Perforated
35842	Beige	3,2	L	Mini Perforated

#### INDICATIONS


Rheumatoid arthritis: Preventive and corrective use for common deviations

#### CONTRAINDICATIONS

This pattern is NOT SUITABLE for the following applications:

- Non-stereotypical deformations
- Very long-term use
- Functional support

#### SIZING



##### SIZING


Measure the total hand width at the MCP joints of the 2nd to 5th finger.

60-70 mm (2 1/5" - 2 3/4") = **SMALL**  
 80-90 mm (3 1/8" - 3 1/2") = **MEDIUM**  
 100-110 mm (4" - 4 3/10") = **LARGE**

#### GENERAL INFORMATION

To take precautionary measures, learn about activation methods , maintenance and waste management read the general Instructions For Use of Orfit Classic.

## APPLICATION METHOD



## PRECAUTIONS BEFORE USE

## ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets, and regulatory information, please visit our website [www.orfit.com](http://www.orfit.com).

### Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries.  
This product is a registered trademark from ORFIT INDUSTRIES BV.



### ORFIT INDUSTRIES BV

Vosveld 9A | B-2110 Wijnegem | Belgium  
T (+32) (0)3 326 20 26  
[welcome@orfit.com](mailto:welcome@orfit.com)

[www.orfit.com](http://www.orfit.com)



Ref. Nr.

Version

Last updated

Revision date

35840-35841-35842 ORFIT CLASSIC

30

1/12/2024

1/12/2026