

Orfit® Classic

Gauntlet immobilization orthosis

GENERAL PRODUCT INFORMATION

This precut is made out of a low-temperature thermoplastic material and can be molded directly on the patient.

PRODUCT RANGE

This precut is available in different sizes.

Article Number	Color	Thickness (mm)	Size	Perforation Type
35830	Beige	2	S	Micro Perforated
35831	Beige	2	M	Micro Perforated
35832	Beige	2	L	Micro Perforated

INDICATIONS


- Acute wrist arthritis
- Wrist sprain
- Radial styloiditis
- Carpal tunnel syndrome
- Rheumatoid arthritis (functional orthosis)

CONTRAINDICATIONS

This pattern is NOT SUITABLE for the following applications:

- Acute fractures of wrist
- Reflex sympathetic dystrophy (RSD)
- Prolonged and aggressive use

SIZING



SIZING

Measure the total hand width at the MCP joints of the 2nd to 5th finger.

60-70 mm ($2\frac{1}{5}'' - 2\frac{3}{4}''$) = **SMALL**


80-90 mm ($3\frac{1}{8}'' - 3\frac{1}{2}''$) = **MEDIUM**

100-110 mm ($4'' - 4\frac{3}{10}''$) = **LARGE**

GENERAL INFORMATION

To take precautionary measures, learn about activation methods, maintenance and waste management read the general Instructions For Use of Orfit Classic.

APPLICATION METHOD



PRECAUTIONS BEFORE USE

ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets, and regulatory information, please visit our website www.orfit.com.

Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries.
This product is a registered trademark from ORFIT INDUSTRIES BV.



ORFIT INDUSTRIES BV

Vosveld 9A | B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

www.orfit.com



Ref. Nr.

Version

Last updated

Revision date

35830-35831-35832 ORFIT CLASSIC

27

1/12/2024

1/12/2026