

Orfit®Classic Intrinsic resting hand orthosis

GENERAL PRODUCT INFORMATION

This precut is made out of a low-temperature thermoplastic material and can be molded directly on the patient.

PRODUCT RANGE

This precut is available in different sizes.

Article Number	Color	Thickness (mm)	Size	Perforation Type
35850	Beige	3,2	S	Non Perforated
35851	Beige	3,2	М	Non Perforated
35852	Beige	3,2	L	Non Perforated


INDICATIONS

- Peripheral nerve palsy
- Brachial plexus
- Flacid hemiplegia
- Reflex sympathetic dystrophy
- Burns
- Quadriplegia
- Volkmann contractures

CONTRAINDICATIONS

This pattern is NOT SUITABLE for the following applications:

- (Acute) Fractures of forearm, wrist and hand
- Prolonged and aggressive use



SIZING


Measure the total hand width at the MCP joints of the 2nd to 5th finger.


60-70 mm (2 $^{1/5}$ " - 2 $^{3/4}$ ") = **SMALL** 80-90 mm (3 $^{1/8}$ " - 3 $^{1/2}$ ") = **MEDIUM** 100-110 mm (4" - 4 $^{3/10}$ ") = **LARGE**


GENERAL INFORMATION

To take precautionary measures, learn about activation methods, maintenance and waste management read the general Instructions For Use of Orfit Classic.

APPLICATION METHOD









PRECAUTIONS BEFORE USE

ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets, and regulatory information, please visit our website www.orfit.com.



Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries. This product is a registered trademark from ORFIT INDUSTRIES BV.



www.orfit.com







 Ref. Nr.
 35850-35851-35852 ORFIT CLASSIC

 Version
 30

 Last updated
 1/12/2024

 Revision date
 1/12/2026