

## Orfit® Classic

### Complete thumb orthosis

#### GENERAL PRODUCT INFORMATION

This precut is made out of a low-temperature thermoplastic material and can be molded directly on the patient.

#### PRODUCT RANGE

This precut is available in different sizes.

Article Number	Color	Thickness (mm)	Size	Perforation Type
35900	Beige	2	S	Micro Perforated
35901	Beige	2	M	Micro Perforated
35902	Beige	2	L	Micro Perforated

#### INDICATIONS


- Basal joint osteo arthritis
- Thumb sprain
- Tendinitis
- Rheumatoid arthritis
- After treatment scaphoide fractures

#### CONTRAINDICATIONS

This pattern is NOT SUITABLE for the following application:

- Acute fractures of wrist or thumb

#### SIZING



##### SIZING

Measure the total hand width at the MCP joints of the 2nd to 5th finger.

60-70 mm (2 <sup>1</sup>/<sub>5</sub>" - 2 <sup>3</sup>/<sub>4</sub>" ) = **SMALL**


80-90 mm (3 <sup>1</sup>/<sub>8</sub>" - 3 <sup>1</sup>/<sub>2</sub>" ) = **MEDIUM**

100-110 mm (4" - 4 <sup>3</sup>/<sub>10</sub>" ) = **LARGE**

## GENERAL INFORMATION

To take precautionary measures, learn about activation methods , maintenance and waste management read the general Instructions For Use of Orfit Classic.

## APPLICATION METHOD



## PRECAUTIONS BEFORE USE

## ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets, and regulatory information, please visit our website [www.orfit.com](http://www.orfit.com).

### Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries.  
This product is a registered trademark from ORFIT INDUSTRIES BV.



**ORFIT INDUSTRIES BV**  
Vosveld 9A | B-2110 Wijnegem | Belgium  
T (+32) (0)3 326 20 26  
[welcome@orfit.com](mailto:welcome@orfit.com)

[www.orfit.com](http://www.orfit.com)



Ref. Nr.	35900-35901-35902 ORFIT CLASSIC
Version	27
Last updated	1/12/2024
Revision date	1/12/2026